



Calcium Hardness /
Кальциевая жесткость
(0 - 500 mg/l)



WATER TESTING
MADE IN GERMANY



Instructions
Инструкции

Reagents

Реактивы




FTS1300

TbsHCH1

Calcium Hardness N° 1 Photometer

TbsPCH2

Calcium Hardness N° 2 Photometer



без моющих
средств!



тщательно
очистить



WATER TESTING
MADE IN GERMANY



CAUTION!

Reagents for water-analysis only!
Do not eat! Keep away from
children! Store cool and dry!
Clean tester properly before
testing!



ВНИМАНИЕ!

Реагенты предназначены
исключительно для анализа воды!
Не употреблять в пищу! Хранить
подальше от детей! Хранить в
прохладном и сухом месте! Перед
тестированием тщательно
очистите тестер!



1




Для тестирования залейте 10 мл воды

10 мл



2




1 X Calcium Hardness № 1 Photometer

3



Растворите таблетку ПОЛНОСТЬЮ!

4



1 X Calcium Hardness № 2 Photometer

5



Растворите таблетку ПОЛНОСТЬЮ!

2 МИН



6



6: Вставьте палочку FTS1300 со шкалой для измерения Calcium Hardness 0 - 500 мг/л

без мощных средств!



тщательно очистить

